



## ClearView C One

**This high-quality video magnifier brings simplicity and high-definition images to everyone who needs magnification. The C-shaped design provides more working space while reading your favorite newspaper, enjoying craftwork, or writing a check - without being restricted in movement.**

If you are looking for an easy-to-use video magnifier, the ClearView C One offers the most commonly used functions: continuous magnification to easily zoom in and out, high-contrast viewing modes for better visibility, large overview for better orientation, and adjustable brightness.

**800-826-4200 • [www.optelec.com](http://www.optelec.com)**





# Optelec ClearView C One

## Key benefits



- Easy-to-use video magnifier with basic controls
- High-definition camera for sharp and vivid images
- Bright 22-inch widescreen display
- Continuous magnification: 2.4x - 27x
- 5 reading modes: photo / black on white / white on black / yellow on blue / yellow on black
- Unique design that fits the interior of your home
- The C-shaped arm gives you more freedom of movement
- More working space to read larger documents more easily
- Space-saving Compact or wide Comfort set-up
- Adjustable brightness, overview function, and position locator
- Wireless Bluetooth LE Remote Control



magnification



speech



daily living aids

**Call or visit today:**

Optelec US  
17757 US Highway 19 N, Suite 560  
Clearwater, FL 33764  
800-826-4200 • [www.optelec.com](http://www.optelec.com)

750631-001 Rev. B

*Simple and effective solutions for the blind and visually impaired*