



NEW!

# Optelec MultiView HD

**Designed for everyday living, the new and improved MultiView HD is a lightweight, high-definition desktop video magnifier that is perfect for near, distance and self-viewing.**

Thanks to the integrated high-definition camera which is lighter and easier to use, the MultiView HD offers better images than ever before. Use it to read, write, view photos and do crafts. Alternatively, tilt and rotate the camera to use it up close like a mirror to apply make-up or view notes on a whiteboard across the room. An auxiliary LED light attached to the camera arm provides additional light when needed.

Simply adjust the magnification, select your preferred viewing mode, or use the snapshot function with an easy to use control pad. The MultiView HD is the all-in-one lightweight solution that will save space on your desk and can easily be moved from room to room using the integrated handle.

**Optelec**



# Optelec MultiView HD

NEW  
DESIGN



## Key benefits

- Improved all-in-one design, perfect for near, distance, and self-viewing
- Magnify text and objects with an easy to use control pad
- High Definition camera with focus lock
- Adjustable document viewing modes for easy reading
  - 2 Photo modes (full color and grayscale)
  - 6 high contrast modes
- Optional reading table for extra support and comfort (sold separately)
- Optional protective travel case with wheels (22" only; sold separately)
- Continuous magnification:
  - 22-inch model: 2X - 75X
  - 24-inch model: 2.4X - 81X
- Dimensions
  - 22-inch model: 21.7" x 17.3" x 7.9"; Weight approx. 17.7lbs
  - 24-inch model: 23.2" x 17.3" x 7.9"; Weight approx. 18.7lbs

*Available through your Authorized Optelec Dealer:*

**Optelec U.S. Inc.**

**Tel: 800.826.4200 | [www.Optelec.com](http://www.Optelec.com)**

**3030 Enterprise Court, Suite C  
Vista, CA 92081**

**Optelec**  
Life is worth enjoying



magnification



speech



daily living aids

**Optelec U.S. Inc.**

Tel: 800.826.4200 - Fax: 800.368.4111 - [info@optelec.com](mailto:info@optelec.com) - [www.optelec.com](http://www.optelec.com)

**Optelec has offices in:**

THE NETHERLANDS – GERMANY – BELGIUM – UNITED KINGDOM – USA – CANADA

*Simple and effective solutions for the blind and visually impaired*