




Optelec ClearView+HD

Regain Your Independence!

Read your morning newspaper, view family photos, enjoy hobbies, or write checks. You can do all this and much more with just “one-touch” of a button.

The ClearView+ is easy-to-use with just one master dial on the table. The dial allows you to enlarge materials to a more comfortable size or change the high contrast viewing modes. Smoothly move the Easy-Glide Reading table in any direction and quickly lock your reading location with the “one-touch” electronic brake. The monitor can be ergonomically positioned, thanks to the ultra flexible telescope arm.





Optelec ClearView⁺

Key benefits



- 5-Directional Ultra Flexible Arm allows the user to adjust the monitor to their preferred position with ease
- Continuous magnification of 1X to 72X
- Full page overview and smooth characters, even at low magnification
- One-Touch electronic lock stops the tray table from moving
- Auto-focus High Definition camera for superior image quality
- Integrated LED lighting for even illumination and elimination of glare
- 22 inch high definition widescreen monitor
- Optional Advanced Feature Pack:
 - Choose from 16 high-contrast color combinations for easier reading
 - Don't lose track again thanks to Position locator and Find function
 - Adjustable window blind and line marker allows to eliminate distractions
 - PC compatible

Available through your Authorized Optelec Dealer:



magnification



speech



braille

Simple and effective solutions for the blind and visually impaired